



# PRINCE OF WALES HOSPITAL

## 威爾斯親王醫院

30 – 32 Ngan Shing Street, Sha Tin, N.T., Hong Kong Tel: (852) 2632 2211 Fax: (852) 2637 8244  
香港新界沙田銀城街三十至三十二號 電話: (852) 2632 2211 傳真: (852) 2637 8244

From: Prof YM Dennis LO  
Hon Chief-of-Service  
PWH Dept of Chemical Pathology

To: All CC, COS, DOM, PWH/SH/BBH/SCH

Tel / Fax No.: 2632 3338 / 2636 5090

Cc: Dr Hong FUNG, HCE(PWH) / CCE(NTEC)  
Dr Susanna LO, SD(CS) / HCE(SH/BBH)  
Dr Michael SUEN, CC Path(NTEC)  
All Duty Biochemists, PWH(CP)  
Ms Candy CHEUNG, CDM(NTEC)  
Ms Karen LAW, DM, PWH(CP)

Date: 18 April 2011

---

### Use of Green-top Heparin Gel Separator Bottle for Iron Profile and Ceruloplasmin Requests

**Implementation Date: 11 May 2011 (6:00 pm)**

For improving operation efficiency, the test mix for 5-mL Green-top Heparin Gel Separator would be expanded to include test requests for Iron Profile (Iron and Total Iron Binding Capacity) as well as Ceruloplasmin. For paediatric patients, 1.3 mL lithium heparin tubes are recommended for the captioned test requests.

The GCRS interface will be updated at 6:00 pm on the implementation date as instructed by Hospital Authority Information Technology Department. Details on sample collection requirements for individual test will be printed on the GCRS specimen labels and can be found in the Test Menu posted at the Prince of Wales Hospital Chemical Pathology i-Hospital homepage ([iCHEMPATH@PWH](mailto:iCHEMPATH@PWH)).

For enquiry, please contact our staff in Specimen Reception at 2632 2363 or 2632 2331.

Thank you for your kind attention.

Sincerely yours,

Prof YM Dennis LO