

\_\_\_\_\_Hospital, NTEC

**(BGA) Internal QC Record**

Serial Number of BGA: \_\_\_\_\_

Expiry Date: \_\_\_\_\_

Opened on: \_\_\_\_\_

To be discarded on: \_\_\_\_\_

1. A new bottle / pack of strip is opened
2. QC testing result is out of the acceptable range.
3. Patient's CBG result is questionable.
4. BGA is dropped.
5. Accuracy of blood ketone result is questioned.
6. Battery of BGA has been changed.

Please change a new form when changing a new BGA / control solution

[illegible]